
new
way
dental
technology

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Gebrauchsinformation / Instruction for use

deutsch | english | français | italiano | español

Stand: 02/2021
Rev. 2

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Fax: +49 (0) 231 - 55 56 - 30

optiprint splint 385
C€ 0044

3D printing resin for additive manufacturing of bite splints (385 nm).

optiprint splint CL plus
C€ 0044

3D printing resin for additive manufacturing of bite splints, drill-guides and orthodontic appliances (385 nm / 405 nm).

optiprint denture
C€ 0044

3D printing resin for additive manufacturing of denture bases (385 nm / 405 nm).

optiprint temp
C€ 0044

3D printing resin for additive manufacturing of temporary crowns and bridges (385 nm / 405 nm).

optiprint guide 385 / 405 / HR
C€

3D printing resin for additive manufacturing of drill-guides (385 nm / 405 nm).

optiprint I-B-T 385
C€

3D printing resin for additive manufacturing of transfer guides for bracket positioning (385 nm).

optiprint tray 385 / HR
C€

3D printing resin for additive manufacturing of customized impression trays (385 nm / 405 nm).

optiprint cast

3D printing resin for additive manufacturing of casting units (385 nm / 405 nm).

optiprint gingiva

3D printing resin for additive manufacturing of gingiva masks (385 nm / 405 nm).

optiprint match	3D printing resin for additive manufacturing of scannable dental form parts (385 nm / 405 nm).
optiprint model	3D printing resin for additive manufacturing of dental models (385 nm / 405 nm).
optiprint model clear 385	3D printing resin for additive manufacturing of demonstration models (385 nm).
optiprint model align	3D printing resin for additive manufacturing of aligner models (385 nm/ 405 nm)

Product description

optiprint are light-curing premium resins based on methacrylate for use in 3D printers (DLP; SLA; CLIP; LCD).

Indication

For additive manufacturing of dental form parts with 3D printing, see details above.

Contraindication

The material should not be used for any other purpose than the additive production of the specified dental parts.

Do not use the polymerized material if allergies exist against any of the ingredients or if contact allergies exist.

Safety instructions

- **optiprint** can irritate eyes and skin (please wear protective goggles).
- Avoid skin contact with non-polymerized material. In rare cases allergic reactions to **optiprint** resins components may occur.
- Use by qualified personnel only.
- Patient target group of the medical devices: adults and adolescents

Risk minimization

- Improper use and non-observance of specifications from these instructions for use will result in quality impairments.
- Keep out of the reach of children!
- During processing, wear protective gloves (nitrile gloves), protective clothing and safety goggles.
- Avoid contact with the eyes. In the event of liquid material making accidental contact with the eyes, immediately rinse eyes thoroughly with plenty of water and consult a physician if necessary.
- Avoid direct contact by the liquid material and ingredients with the skin before rehardening.
- Biocompatibility is only guaranteed with proper use (polymerization in a protective atmosphere).
- Please observe the safety data sheet.

Directions for Use

To achieve optimum characteristics, please proceed as follows:

Before and during the printing process:

- **optiprint** bottles should be swiveled several times (not shaken) before use.
- Pour **optiprint** carefully into the specified 3D printer tray.
- Allow the material to rest in the tray until all bubbles in the material have disappeared.
- Observe a processing temperature of 23°C - 40°C.
- The 3D printer should be in working order, clean and adequately filled with **optiprint**.

After printing

- Recommended drip off time after booting up the platform: appr. 10 minutes
- Immediate post-processing is recommended afterwards. If this cannot be guaranteed, leave the produced objects suspended in the printer.
- Clean the objects with suitable cleaning solution (2x 3 min in isopropanol - pre-cleaning/post-cleaning - in the ultrasonic basin)
- The dental form parts produced with **optiprint** can be polished and repaired in the usual manner.
- **optiprint** can be filtered thanks to its low viscosity. Recommended filter diameter: <100µm.
- Fill the **optiprint** material from the tray back into the bottle within 12 hours after use.

Post-polymerization

Product	optiprint: temp, denture, guide 385 / 405 / HR, I-B-T 385, tray 385 / HR, gingiva, model, model clear 385, model align, match, splint 385 / CL Plus	optiprint: cast
Light output in the Otoflash G171	2x 2000 flashes (turn components around after 2000 flashes)	1x 750 flashes (Object on the model)

General information

- Ensure a clean job. Impurities on the machine may cause errors on the printed object and damage the tray.
- Comply with machine and material parameters.
- Minimum material thickness for splints with **optiprint** splint and **optiprint** splint CL plus: 1.5 mm to ensure increased resistance to breakage.
- Optimal layer thickness for all **optiprint** Class IIa products: 50 µm
- Maximum full hardening depth during post-exposure: For solid objects and exposure on both sides, the material thickness may be up to 7 mm (at a full hardening depth of 3.5 mm).
- Pay attention to the operating instructions of the 3D printer
- Observe the parameters of the printer software

Composition

Mixture of acrylate and methacrylate resins, photoinitiators (phosphine oxide), additives, dyes.

Storage and shelf life

- Store the close units at room temperature (18°C - 28°C) in a dark and dry place. Even slight exposure to light may trigger polymerization.
- Remark: When handling **optiprint**, only expose it to daylight/ambient light for as short a period as possible, because if the material is exposed to light for too long, this will result in inadvertent curing.
- Always keep the bottle tightly sealed. After using, immediately reseal carefully.
- Do not use after the expiry date!

Information

dentona AG is certified in accordance with DIN EN ISO 13485 and operates an elaborate quality assurance system. All dentona products are thus free of material and manufacturing defects and always have flawless quality. Our user recommendations are based on a careful determination of practice-related values from our dental laboratory. Achieving the guaranteed quality and usability of the product requires strict compliance with all process steps specified in the product files in a two-tier quality control. The user is responsible for operation and proper use of the material. The medical products are intended for use as custom-made products by dental technicians in accordance with EU regulations.

Report all serious incidents (death, serious deterioration of health, serious risk to public health) that have occurred in connection with the medical device to the manufacturer and the competent authority of the member state.

If defects appear in the material within the warranty period, the user only has a right to replacement of the material. dentona AG is liable for neither losses nor damage caused by this material, irrespective of whether this damage is direct or indirect, and especially collateral or consequential damage, regardless of the legal basis. According to statutory regulations, dentona AG is liable for direct damage to the material based on a deliberate act or gross negligence by its legal representatives or executive staff and for personal damage. Any liability for the material and damage as a consequence of using it are excluded if the user has not observed the specified process steps.

Disposal

Dispose according to official regulations. Must not be disposed of together with household waste. Do not empty into drains.

Symbols

	See instruction for use (also on www.dentona.de)		temperature limitation
	Use by/Expiry date		Lot number
	Item number		Keep away from sunlight
	Manufacturer		CE mark (with or without number of notified body)
	Attention health hazard		Attention environmental hazard
	Attention corrosive		Attention serious health hazard
	Medical device		

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